







# XploRA Series

Simply Better Raman



Chemical Analysis & Imaging  
for Research & Industry

# Why Raman ?

## Advantages


- Non-destructive
- Non-contact/*In-situ* sampling
- Reduced sample preparation
- Water/aqueous phase sampling
- Organic/inorganic molecules
- Amorphous/crystalline

Raman is an ideal technique for research and industry offering high quality data, reliability, versatility and improved value for money over other analytical techniques. Benefits not only include the range of samples that are suitable for analysis, but also the information content that is provided.

- Chemical identification
- Quality testing
- Process/product troubleshooting
- Contamination and inclusion analysis
- Raw materials inspection

## Principle

Interaction of laser light with a sample results in a Raman spectrum - a detailed chemical fingerprint. Combined with an optical microscope, this provides sample identification and chemical imaging on a microscopic scale.




The sensitivity, spectroscopic and imaging performance of the XploRA product line

enables the broadest range of sample analysis.


## Pharmaceuticals




## Nano-materials



## Semiconductors



## Polymers




## Geology




## Chemicals



## Art and Museum



## Forensics



# XploRA Series



## XploRA PLUS: Research

➤ Raman imaging has never been so fast!

The XploRA PLUS incorporates unique and powerful functions in a reliable, high performance system, ideally suited to the research and analytical lab.

It is **fully confocal**, not compromising image quality, spatial or depth resolution. The **SWIFT™ Fast Raman images** are the fastest fully confocal Raman images available, typically 10x faster than conventional Raman imaging.

The simplicity and power of the XploRA PLUS is unmatched with an enhanced range of options such as multiple laser wavelengths, EMCCD detection, Raman polarisation and even AFM Raman TERS combination.

It is the best platform for multi-sample and multi-user environments.

- Fastest confocal imaging
- Automated laser wavelength switching with just a single mouse click
- Large range of options and accessories



## XploRA ONE: Analytical

➤ Raman analysis has never been so easy!

The XploRA ONE offers new capabilities to the industrial and analytical user, providing the highest performance Raman, in a **cost effective and robust instrument package**.

It is ideal for routine analytical, research and quality testing applications.

- OneClick operation
- Auto-calibration in OneCheck
- Regulatory compliance: 21 CFR11



## XploRA INV: Life Science

➤ Hybrid biological imaging and analysis made easy!

The XploRA INV Raman microscope is the only **truly analytical inverted Raman microscope**. Configured for high sensitivity bio-Raman spectroscopy. It offers TRUE confocal performance with low maintenance and dedicated software for ease of operation.

The uniquely integrated system design ensures stability, optimizing the imaging workflow. The integrated inverted microscope enables multi-modal analysis, such as **fluorescence, Raman, laser tweezing and even TERS analysis** to be conducted upon the same instrument and at the same sample position.


- Inverted life-science microscope
- Fast and simple sample imaging
- Multi-modal analysis



# The XploRA

Innovation provides improved productivity and extended reliability.

- Superior CCD and EMCCD for Class leading S:N
- Easier and faster analysis
- Requiring far less laser power on the sample, preserving sample integrity
- Ability to detect thin films, small particles, and dilute solutions




Silicon (Si) 4th order sensitivity

## Improved Detection and Sensitivity

## Raman Imaging just got *Faster & Easier* with SWIFT™ and SWIFTXS


- Suitable for all laser wavelengths
- Scalable to large area and micron scale imaging for maximum image detail
- Means faster and more reliable Raman imaging at the click of a button
- Fully confocal for improved image detail

Combined with an EMCCD, **SWIFTXS** pushes the speed of confocal Raman imaging to the next level. Accelerate your chemical imaging and generate high definition Raman images in minute timescale. **Supercharge your Raman Imaging!**



Fast SWIFT™ image of semiconductor defect, obtained in less than 2 minutes

Particulate contamination Raman SWIFT™ Image



40 200 spectra acquired in less than 50 seconds  
Structured semiconductor device, Raman SWIFTXS image




# At a Glance

Technical design offers automation and class leading ease of operation.



## Full Optical Microscopy To See Your Sample




Extended polarized light microscope image

- All modes of microscopy, transmission and reflection illumination
- Options for DIC, phase, epifluorescence, dark field and polarized light microscopy
- Automated Extended Video Montage
- Range of options such as Autofocus, and ParticleFinder (auto-location) will never limit the scope of the optical microscopy

## Maximum Detail, Resolution and Range

- Optimized range and resolution in one-shot for all lasers
- Full range optics enables detail over the complete Raman spectrum
- The high resolution ability of the XploRA PLUS and INV, provide the highest spectroscopic detail while optimizing the Raman sensitivity



Raman spectra of pharmaceuticals and organics can often show subtle information in the 2800-3400  $\text{cm}^{-1}$  region

# Raman Data in Seconds: For Industrial Quality Control, Research and Analytical Testing

The XploRA Series offers full spectral analysis in OneClick Operation.



Simple operation, speed, low maintenance and push button results enables you to optimize productivity and efficiency.

The XploRA Series is driven by LabSpec's intuitive user interface enabling logical workflow through an experiment.

- HORIBA OneClick Raman operation
- Autocalibration
- Enhanced multi-page analysis reports
- Augmented help and troubleshooting
- Sample methods for routine repetition
- User login accounts for system security
- 3D volume and topography imaging
- Enhanced on-the-fly image generation
- Support of multi-screen PC environments
- Extended microscope image correction
- Extensive spectral database/libraries
- Easier chemometric processing
- ParticleFinder for automated particle location



## Easily Expand your Analytical Technique

Benefits of XploRA Raman Microscopy in Research and Industrial Analysis		Compared to
<ul style="list-style-type: none"> <li>• No sample preparation</li> <li>• Full optical microscope</li> <li>• Water based samples</li> <li>• Sub-micron scale analysis</li> </ul>	<ul style="list-style-type: none"> <li>• Non-contact analysis</li> <li>• Multi-layer analysis</li> <li>• Inorganics</li> <li>• Polymer backbone characterization</li> </ul>	FTIR
<ul style="list-style-type: none"> <li>• Non-destructive analysis</li> <li>• No sample preparation</li> <li>• Solids/surface analysis</li> </ul>	<ul style="list-style-type: none"> <li>• Molecular/crystal structure</li> <li>• Fast analysis times typically &gt; 2 seconds</li> <li>• Low maintenance</li> </ul>	Mass Spectrometry
<ul style="list-style-type: none"> <li>• Crystalline and amorphous materials</li> <li>• Single particle analysis</li> </ul>	<ul style="list-style-type: none"> <li>• In-situ environmental stages (heating/cooling/relative humidity)</li> <li>• Small benchfootprint</li> </ul>	XRD
<ul style="list-style-type: none"> <li>• Chemical information and chemical images</li> </ul>		Optical Microscopy
<ul style="list-style-type: none"> <li>• No sample preparation</li> <li>• Environmental conditions: no vacuum, heating/cooling/relative humidity controlled stages</li> <li>• Fast analysis times and sample throughput</li> </ul>	<ul style="list-style-type: none"> <li>• Fast start-up and ready to analyze time in less than 10 minutes from off</li> <li>• Multi-layer samples</li> <li>• Small bench footprint</li> </ul>	SEM

# Simplified Workflow

With  **LabSpec 6**  
Spectroscopy suite



Auto-Calibration  
OneCheck

Input Sample  
No preparation

Focus  
Microscope

Easy  
Navigation

Raman  
Acquisition

Obtain Result  
Match to Library

Report Result

## Simplified Calibration and Validation

The patented calibration objective tool uses a certified reference material and ASTM method for system validation and ensures fast, easy calibration and validation of the instrument.




OneCheck and the system is ready to run.

Impressively fast start-up time, less than 10 minutes from cold, removes the need for any lengthy start-up procedures, adjustment or the need for the system to be continually powered on. Ideal for laboratory efficiency, running costs and the environment!

## EasyNav™



With the groundbreaking EasyNav™ package, it is fast and easy to navigate in-focus, in real-time, to identify the region of interest and obtain sharp, clear Raman chemical images, thanks to three new revolutionary applications: NavMap™, NavSharp™ and ViewSharp™.



NavMap™ View

## OneClick Raman Acquisition



OneClick Raman acquisition optimizes acquisition parameters and signal processing in OneClick, including baseline corrections, fluorescence rejection and noise reduction.

# XploRA Series Specifications

	Industrial	Research	Bio/Life science	Nano Raman
	XploRA™ One	XploRA™ PLUS	XploRA™ INV	XploRA™ AFM Raman
Faster Raman SWIFT™ Imaging / SWIFT XS (with EMCCD)	YES with XY stage	YES with XY stage	YES with XY stage Optional DuoScan™ imaging	YES Optional XY stage & DuoScan™ depending upon version
Confocal Imaging	1 µm XY	0.5 µm XY	0.5 µm XY	0.5 µm XY 10 nm* with TERS
Routine operation Automation	OneClick Auto/Raw	OneClick Auto	Methods and scripts	Integrated AFM Raman software
Full Microscope	Upright	Upright	Inverted	Upright and/or inverted
Resolution	Standard	Standard + High	Standard + High	Standard + High
Multi-laser Options	Single laser 532 and 785 nm	532, 638, 785 nm others on request	532, 638, 785 nm others on request	532, 638, 785 nm others on request

\* Requires HORIBA TERS tips

## XploRA Nano

- AFM Raman module
- XploRA PLUS and INV versions
- TERS Ready: 10 nm resolution\*
- Multi-sampling geometry: upright-inverted-side axis
- SWIFT™ Nano-Raman images
- High performance AFM functionality
- Integrated control and construction
- Stability and reliability

The XploRA PLUS and INV can add the NanoRaman (TERS) extension to probe nanometer structures and single molecules in a single compact, high performance system.



$\lambda = (400 \text{ nm} - 800 \text{ nm})$   $P \leq 150 \text{ mW}$   
**VISIBLE AND/OR INVISIBLE LASER RADIATION**  
**AVOID EXPOSURE TO BEAM**  
**CLASS 3B LASER PRODUCT**



Find out more at [www.horiba.com/xplora](http://www.horiba.com/xplora)

[www.horiba.com/scientific](http://www.horiba.com/scientific)  
[info.sci@horiba.com](mailto:info.sci@horiba.com)

**HORIBA**  
 Scientific

**France:** HORIBA Jobin Yvon S.A.S., 16-18 rue du Canal, 91165 Longjumeau cedex - Tel. +33 (0)1 69 74 72 00 - Fax. +33 (0)1 69 09 07 21 - Email: [info-sci.fr@horiba.com](mailto:info-sci.fr@horiba.com)  
**USA:** HORIBA Instruments Inc., 20 Knightsbridge Road, Piscataway, NJ 08854 - Tel. +1 732 494 8660 - Fax. +1 732 549 5125 - Email: [info-sci.us@horiba.com](mailto:info-sci.us@horiba.com)  
**Japan:** HORIBA Ltd., 2 Miyahogigashi, Kisshoin, Minami-ku, Kyoto 601-8510 - Tel. +81(75)313-8121 - Fax. +81(75)321-8312 - Email: [info@horiba.co.jp](mailto:info@horiba.co.jp)  
**Germany:** HORIBA Jobin Yvon GmbH, Neuhofstrasse 9, 64625 Bensheim - Tel. +49 (0) 6251 8475 0 - Fax. +49 (0) 6251 8475 20 - Email: [info-sci.de@horiba.com](mailto:info-sci.de@horiba.com)  
**Italy:** HORIBA Jobin Yvon Srl., Viale Luca Gaurico 209 - 00143 Roma - Tel. +39 06 51 59 22 1 - Fax. ++39 06 51 96 43 34 - Email: [info-sci.it@horiba.com](mailto:info-sci.it@horiba.com)  
**UK:** HORIBA UK Ltd., Kyoto Close, Moulton Park, Northampton, NN3 6FL - Tel. +44 (0)1604 542 500 - Fax. 44 (0)1604 542 699 - Email: [info-sci.uk@horiba.com](mailto:info-sci.uk@horiba.com)  
**China:** HORIBA (China) Trading Co. Ltd., Unit D 1F, Bldg A, Srynnex International Park, No. 1068 West Tianshan Road, Shanghai 200335 - Tel. +86 (0)21 6289 6060 - Fax. +86 (0)21 6289 5553  
 Email: [info-sci.cn@horiba.com](mailto:info-sci.cn@horiba.com)  
**Singapore:** HORIBA Instruments (Singapore) Pte Ltd. 3 Changi Business Park Vista #01-01 AkzoNobel House 486051 Singapore - Tel. +65 (6) 745-8300 - Fax. +65 (0)6 745 8155  
**Taiwan:** HORIBA Taiwan Inc., 8F-8, No.38, Talyuan St., Zhubei City, Hsinchu County 30265 - Tel. +886 3 5600606 - Fax. +886 3 5600550 - Email: [info-sci.tw@horiba.com](mailto:info-sci.tw@horiba.com)  
**India:** HORIBA India Private Limited Bangalore Office No.55, 12th Main, Behind BDA Complex, 6th sector, HSR Layout, Bangalore South, 560102 Bangalore Tel. +91 (80)4127 3637  
**Brazil:** HORIBA Instruments Brasil Ltda., Rua Presbítero Plínio Alves de Souza, 645, Loteamento Polo Multivias Bairro Jardim Ermida II, Jundiá São Paulo CEP 13.212-181  
 Tel. +55 (0)11 2923 5400 Fax. +55 (0)11 2923 5490 - Email: [infocientifica.br@horiba.com](mailto:infocientifica.br@horiba.com)  
**Other:** Tel. +33 (0)1 69 74 72 00 - Email: [info.sci@horiba.com](mailto:info.sci@horiba.com)