




## Applications

- Semiconductor characterization
- Photovoltaic Solar Cell Testing
- UV Exposure Testing
- Sunscreen Testing
- Cosmetics Testing
- Environmental Testing
- Electrochemical

## Features

- Economical Design
- Up to Class AAA Specification
- Touchscreen Power Supply
- Turn Key Operation
- Collimated Systems Available
- Manual Shutter Included
- Electronic Shutter Optional
- Multiple Optional Accessories
- Lamp Life Timer
- Air Mass AM0/AM1.5 Filters

# SF SOLAR SIMULATOR

## Small Area Collimated Lens Based

## Class AAA and ABA

# SF Solar Simulator


## OVERVIEW

Sciencetech's SF Solar Simulators (Steady State) are low cost lens based systems designed for researchers who do not require a large field of illumination. SF series solar simulators produce 1 Sun and are available in Class A or B uniformity.

The beam can be projected horizontally (standard) or vertically with the use of a downward-facing stand.

Sciencetech SF series solar simulators produce a highly collimated output and are an ideal choice for space based research or systems needed high levels of collimation.

Sciencetech SF type Solar Simulators include an arc lamp housing, one xenon arc lamp, touchscreen power supply with igniter, filter holder, and testing report.





Vertical output with downward facing stand

Horizontal output


## Standards for Class AAA Specifications

Sciencetech's solar simulator specifications listed are according to ASTM E927 standards, unless otherwise stated. Please contact us if you are interested in matching IEC 60904-9(2007), JISC 8912-1998, or other standards. We can accommodate testing to match several standards.


**Class A. Spectral Match.** Solar simulator spectrum meet with ASTM AM1.5G and AM0 solar spectrum.



**Class A. Non Uniformity of SF300A over 1" Diameter** less than 2%



**Class A. Temporal Instability of Irradiance.** Less than 2%.



# SF Solar Simulator

## SPECIFICATIONS

### SF Series Models

| Model                                               | SF                                                        |            |
|-----------------------------------------------------|-----------------------------------------------------------|------------|
| Part Number                                         | 300                                                       | 150        |
| Solar Simulator Classification                      | AAA to UBA                                                | ABA to UBA |
| Spectral Range (nm)                                 | 250-2000                                                  |            |
| Spectral Match Classification <sup>1</sup>          | Up to A                                                   |            |
| Spatial non-uniformity                              | < 2%                                                      | < 5%       |
| Non-Uniformity Class                                | A                                                         | B          |
| Temporal Stability Classification                   | A                                                         |            |
| Target Diameter (mm)                                | 25 to 50                                                  | 50         |
| Working Distance (mm)                               | 100-130                                                   |            |
| Collimation                                         | 1.0 degree half angle <sup>2</sup>                        |            |
| Power Level at Target (AM1.5G Standard — 100mW/cm²) | 1 Sun                                                     |            |
| Center Beam Line Height (mm)                        | 137                                                       |            |
| Lamp Power (W)                                      | 300                                                       | 150        |
| Power Supply Model                                  | 601-300                                                   | 601-150    |
| Dimensions (LxWxH) (mm)                             | 222 x 183 x 318                                           |            |
| Weight (kg) (with power supply)                     | 6 (12)                                                    |            |
| Power Supply Input                                  | 110-240V, 50Hz/60Hz, 450W                                 |            |
| Operating Current (A)                               | 5-20                                                      |            |
| Stability / Ripple / Regulation                     | 0.05% / < 1% / 0.02% current variation for 5V line change |            |

1. Filter required for spectral match classification
2. >70% power within specification



[Browse Solar Filters](#)

Sciencetech's low cost line of SF solar simulators include a filter box which can hold a range of filters in Sciencetech's standard FT style filter holder. The most popular options are AM filters; however, a range of other filter options are available.

**AM0 Filter** reproduce extra-terrestrial solar spectrum, used for space applications. **AM0 spectral match only.**

**AM1.5 Global** simulates the global total radiation solar on the ground when the sun is at 48.2° zenith angle. It includes both direct light from the sun and the diffuse light that is scattered by the atmosphere.

**AM1.5 Direct.** Reproduce the direct radiation spectrum on the ground at 48.2° zenith angle.

### Sciencetech Filters

| Model    | Description                                                 |
|----------|-------------------------------------------------------------|
| 160-8090 | (AMF-AM1.5G) Air Mass AM1.5G Filter for SF Series           |
| 160-8091 | (AMF-AM1.5D) Air Mass AM1.5D Filter for SF Series           |
| 160-8092 | (AMF-AM0) Air Mass AM0 Filter for SF Series                 |
| 100-8048 | (WF-1Q) Compact IR water Filter, 1.75 " with Quartz Windows |



[Browse all Filtering Options](#)


# SF Solar Simulator CONFIGURATION

## Power Supply

Sciencetech's 601 – series power supplies are the included power supplies for use with Sciencetech's SF series lamp houses.

**Standard features** included with Sciencetech's 601– series power supplies:


- Touchscreen interface
- Shutter and exposure control (if electronic shutter is supplied\*)
- Single connection for lamp power, cooling, and communication
- Lamp starts and timer log
- Fan cooling safety interlock
- RS232 software GUI included shown below



## Optional Upgrades

To be added to sales order as optional upgrades:

- Optical feedback
- Auto lamp starting




Contact a Sciencetech Technical Sales Representative  
to discuss your custom requirements!

# SF Solar Simulator

## ACCESSORIES


Sciencetech manufactures modular spectroscopy and solar simulation equipment. The SF type simulators are based on Sciencetech's compact LH series lamp house; due to this modular design philosophy, there are a number of available options for SF style solar simulators from Sciencetech's catalog of instrument accessories.



### Stand (101-8015)

---

Downward facing stand for LH series lamp houses.



### Automated Shutter 2" (127-9005)

---

Computer controlled shutter for LH series lamp houses (works with SF series solar simulators).



### Solar Reference Cell (125-9060)

---

This solar reference cell consists of a 20 × 20 mm monocrystalline silicon photovoltaic cell encased in a 92 × 70 × 16 mm metal enclosure with a protective quartz window and a temperature sensor. Options are available for calibration by Sciencetech, or by ISO/IEC 17025 accredited lab. For ISO/IEC 17025 accredited lab calibrated cells, calibration is traceable to NIST and NREL.



### IV tester (175-9120)

---

Current Voltage Measurement system (IV Tester) for Continuous Solar Simulators.



### Dark safety glasses (720-0159)


---

UV Dark safety glasses.

# SF Solar Simulator

## DIMENSIONS

Dimensions are in inches and [mm].



# SF Solar Simulator

## ORDERING INFORMATION

| Model           | Part Number  | Description                                                                |
|-----------------|--------------|----------------------------------------------------------------------------|
| SF150-1.5G-ABA  | 160-9050-C05 | 150 W Arc Lamp, AM1.5G Filter, Class ABA, 1" Diameter Target               |
| SF150-1.5D-ABA  | 160-9050-C06 | 150 W Arc Lamp, AM1.5D Filter, Class ABA, 1" Diameter Target               |
| SF150-0.00-ABA  | 160-9050-C07 | 150 W Arc Lamp, AM0 Filter, Class ABA, 1" Diameter Target                  |
| SF150-X.XX-UBA  | 160-9050-C08 | 150 W Arc Lamp, No Filter, Class UBA, 1" Diameter Target                   |
| SF300-1.5G-AAA  | 160-9050-C09 | 300 W Arc Lamp, AM1.5G Filter, Class AAA, 1" Diameter Target               |
| SF300-1.5D-AAA  | 160-9050-C10 | 300 W Arc Lamp, AM1.5D Filter, Class AAA, 1" Diameter Target               |
| SF300-0.00-AAA  | 160-9050-C11 | 300 W Arc Lamp, AM0 Filter, Class AAA, 1" Diameter Target                  |
| SF300-X.XX-UAA  | 160-9050-C12 | 300 W Arc Lamp, No Filter, Class UAA, 1" Diameter Target                   |
| SF300-1.5G-ABA  | 160-9050-C13 | 300 W Arc Lamp, AM1.5G Filter, Class ABA, 2" Diameter Target               |
| SF300-1.5D-ABA  | 160-9050-C14 | 300 W Arc Lamp, AM1.5D Filter, Class ABA, 2" Diameter Target               |
| SF300-0.00-ABA  | 160-9050-C15 | 300 W Arc Lamp, AM0 Filter, Class ABA, 2" Diameter Target                  |
| SF300-X.XX-UBA  | 160-9050-C16 | 300 W Arc Lamp, No Filter, Class ABA, 2" Diameter Target                   |
| HPF-LP-400-FT-3 | 640-9021     | Long Pass Filter, 400nm cut-on                                             |
| HPF-LP-420-FT-3 | 640-9020     | Long Pass Filter, 420nm cut-on                                             |
| HPF-LP-NIR-FT-3 | 640-9015     | Long Pass Filter, 700nm cut-on                                             |
| HPF-LP-UVA-FT-3 | 640-9017     | Long Pass Filter, 320nm cut-on                                             |
| HPF-LP-VIS-FT-3 | 640-9014     | Long Pass Filter, 400nm cut-on, for mid-range power                        |
| HPF-BB-NIR-FT-3 | 640-9012     | Hot Mirror for UV Applications, 700-1100nm block                           |
| HPF-BB-UVB-FT-3 | 640-9018     | Bandblock Filter, 300-420nm block                                          |
| HPF-BB-VIS-FT-3 | 640-9016     | Bandblock Filter, 400-700nm block                                          |
| HPF-BP-VIS-FT-3 | 640-9019     | Bandpass Filter, Visible, Blocks UV/NIR                                    |
| MF-FT-3         | 640-9050     | Mounted Mesh Neutral Density Filter, Various Transmissions                 |
| Custom Filters  | Custom       | Custom filters—please inquire.                                             |
| LH-DFS          | 101-8015     | Downward Facing Stand for LH series and SF Solar Simulators                |
| SOL-METER       | 125-9011     | Calibrated Sun Meter                                                       |
| SOL-REF-Q       | 125-9060     | Calibrated Reference Cell, Sciencetech or external lab calibration options |
| SCI-MO          | 125-9040     | Monitor Detector Assembly                                                  |
| MULTIMETER      | 125-9030     | Multi-ranging Multimeter and Test Leads                                    |
| SSIVT-C         | 175-9120     | Current-Voltage Measurement System (IV Tester) for Continuous Solar        |
| SCP-4G          | 165-8214     | Probe Station, 4 Probes, Gold Spring Probes                                |
| SCI-SCC3-TE     | 165-8202     | Solar Cell Chuck, TE-cooled, Vacuum-ready, 9 cm × 9 cm (3.5" × 3.5")       |
| SCI-SCC3-L-B    | 165-8221     | Solar Cell Chuck, Liquid-cooled, Rear Contact, 9 cm × 9 cm (3.5" × 3.5")   |