

Laser Safety Doggles® Kit



Articlenumber: DKTP5K015000
GTIN: 4050369018266
Unit: 1 Stück
Weight incl. packaging 0,23 kg
Weight excl. packaging: 0,12 kg

Highlights

The Laser Safety Doggles® DKTP5K015000 with laser safety filter P5K01 protect the eyes of dogs, cats, kangaroos and other little animals dependably from NIR-laser radiation (800-1100nm, OD5) during laser procedures. The elastic headband ensures a secure and dependable fit during the laser therapy. The Laser Safety Doggles® is available as a Kit. This includes three different sizes.

The Laser Safety Doggles® are recommended by the [Aussie Mobile Vet](#) in Australia which uses them daily during laser therapy sessions with animals.

FILTER:	P5K01
FILTER COLOUR:	Dark green
FILTER MATERIAL:	Plastic
FILTER TECHNOLOGY:	Absorption filter
FILTER THICKNESS:	ca. 2mm
FRAME:	Doggle
FRAME TYPE:	Pets goggle
PROTECTION RANGE:	near infrared
VLT (APPROX.):	6%
VISUAL BRIGHTNESS:	Limited
COLOUR RECOGNITION:	Limited